

Aesthetic Denture Wax



Provisional wax aesthetics

Provisional wax aesthetics

The new Aesthetic Denture Wax and Aesthetic Color Wax-Set make it possible to give the patient a true-to-life picture of the finished denture even during the wax try-in.



Aesthetic Denture Wax



Everyday wax

Aesthetic Denture Wax SOFT is noted for its special consistency and specially formulated for everyday use in shade 34 for natural aesthetics and to match the Candolor acrylic.

Advantages

- Shaded in shade 34
- Special malleability
- Ideal flow properties
- Good modeling characteristics
- Easily reduced



Shade 34

Modeling wax

Aesthetic Denture Wax MEDIUM and HARD are special modeling waxes with very low shrinkage, making them ideal for carving, in shade 34 for a natural appearance and to match the Candolor acrylic.

Advantages

- Shaded in shade 34
- No distortion even during a long try-in
- High degree of hardness and solidification point
- Two degrees of hardness (medium and hard)
- Minimal shrinkage
- Optimum modeling properties

Aesthetic Color Wax-set

Characterization wax

In order to make individual characterization of the wax model simpler and faster for the dental technician, the specialists at Candolor have developed three ready-mixed modeling waxes in shades to match the Aesthetic Color Set Easy. For maximum naturalness, shade 53 should be used for the marginal gingiva, 55 for the fixed gingiva and 57 for the alveolar mucosa.



Benefits

- Color coordinated with the Color Set Easy acrylics
- Very hard wax with low shrinkage ideal for carving
- Gives the patient the option of choosing a characterized denture for the final product during the wax try-in
- High gloss after polishing

Candolor AG

CH-8602 Wangen/ZH, Pünten 4, Postfach 89
Tel. +41 (0)44 805 90 00, Fax +41 (0)44 805 90 90
www.candolor.com, candolor@candolor.ch

Subsidiary: Candolor Dental GmbH

D-78239 Rielasingen-Worblingen, Am Riederngraben 6
Tel. +49 (0)7731 79 783-0, Fax +49 (0)7731 28 917
www.candolor.de, info@candolor.de

Subsidiary: Candolor USA Inc.

7462 N. Figueroa St., Suite 104, Los Angeles CA 90041
Phone +1 (323) 254-1430, Phone +1 800 436-3827
Fax +1 (323) 254-5146
www.candolor.us, info@candolor.us